



**Vivex DBM 100™ is a demineralized bone matrix paste derived from 100% allograft bone.**

Vivex DBM 100 exhibits the essential elements of bone formation. It does not contain any extrinsic carriers and is entirely derived from 100% allograft bone. The product is provided in a ready to use syringe.

Indicated for use as a bone void filler for the treatment of surgically created or traumatic skeletal defects that are not intrinsic to the stability of the bony structure.

Vivex DBM 100 may be used in surgical procedures including:

- Spine
- Orthopedics
- Reconstruction
- Dental
- Neuro
- Trauma
- Foot & Ankle

**Features and Benefits:**

- Osteoinductive - each lot is assayed
- Osteoconductive
- 100% bone - no additional carrier
- Bio-compatible
- Quality control tested for cohesiveness, malleability, and flowability

**Ordering information:**

ITEM NUMBER	PRODUCT DESCRIPTION
DBM001	1cc
DBM002	2cc
DBM005	5cc
DBM010	10cc



To place an order, please contact: [customer care@vivex.com](mailto:customer care@vivex.com)  
1951 NW 7th Ave, Suite 200, Miami, FL 33136 | 888-684-7783 | [www.vivex.com](http://www.vivex.com)

Trademarks ™ and Registered Trademarks ® of Vivex Biomedical, Inc. Copyright © 2017 Vivex Biomedical, Inc. All rights reserved.

MKG-30012  
Rev. 04